

## **PREPARATION OF SRM 967 CREATININE IN HUMAN SERUM**

**Prepared at Solomon Park Research Institute, Kirkland, Washington**

**Human serum with the following characteristics:**

**Master Pool**

**Comprised of units drawn from postmenopausal female donors to yield pool with creatinine of 0.8 mg/dL or lower.**

**Collected and handled following NCCLS C37-A**

**All units tested and found negative for viral markers as specified by FDA**

**Pooled and split into two sub-pools**

**PREPARATION OF SRM 967 *continued***

**For Low level, 1200 1 mL aliquots prepared in 2 mL amber vials**

**For High level, second sub-pool spiked with crystalline creatinine to bring the concentration to 4.0 mg/dL  $\pm$  0.2mg/dL.**

**1200 1 mL aliquots prepared in 2 mL amber vials**

## **LC/MS METHOD FOR SERUM CREATININE**

**Based on published method: [P. Stokes and G. O'Connor, J.,  
Chromatog. B 794, 125-136 (2003)]**

### **Method in brief:**

**Spike serum with creatinine-d<sub>3</sub>**

**Precipitate proteins with cold ethanol**

**Centrifuge**

**Decant supernatant and dry under nitrogen**

**Reconstitute in water and filter**

**Dilute with 10 mM ammonium acetate**

**LC/MS**

**Phenomenex LUNA C-18**

**Gradient: 10 mM ammonium acetate for 7 min,**

**Then acetonitrile:10 mM ammonium acetate (80:20) and hold for 13 min**

**Electrospray ionization – positive mode monitoring (M+H)<sup>+</sup> at 114/117**

**Calibration from standard curve using SRM 914a and creatinine- d<sub>3</sub>**

## MEASUREMENT OF SRM 967 USING LC/MS

**3 Sets, 3 Vials of Each Level per Set, 2 Aliquots/vial  
(mg/dL)**

Set	Box	Vial	Sample	Concentration	Set Mean	Set St. Dev.	Set CV, %
SRM 967 Level I							
1	3	39	1	0.749	0.747	0.001	0.180
1	3	39	2	0.748			
1	13	7	3	0.746			
1	13	7	4	0.748			
1	18	49	5	0.746			
1	18	49	6	0.746			
2	4	12	15	0.747	0.747	0.002	0.251
2	4	12	16	0.750			
2	14	22	17a	0.749			
2	14	22	18a	0.745			
2	15	39	19	0.747			
2	15	39	20	0.745			
3	6	13	29	0.749	0.748	0.002	0.315
3	6	13	30	0.750			
3	11	44	31	0.746			
3	11	44	32	0.750			
3	22	28	33	0.746			
3	22	28	34	0.745			
Mean				<b>0.747</b>			
St. Dev.				<b>0.002</b>			
CV, %				<b>0.243</b>			

### **SUMMARY RESULTS FOR LC/MS**

<b>Level</b>	<b>Mean</b>	<b>CV%</b>	<b>Certified Value</b>
<b>967 level I</b>	<b>0.747</b>	<b>0.24</b>	
<b>967 level II</b>	<b>3.918</b>	<b>0.24</b>	
<b>909b level I</b>	<b>0.632</b>	<b>0.67</b>	<b>0.6355 ± 0.0062</b>
<b>909b level II</b>	<b>5.275</b>	<b>0.23</b>	<b>5.287 ± 0.060</b>

## **GC/MS MEASUREMENTS OF SRM 967**

**NIST ID-GC/MS “Definitive Method”**

**Measurements are starting in August**

**New staff member is responsible for setting up method**

## **PROPOSED COMMUTABILITY STUDY**

### **Materials**

**20 fresh frozen patient specimens**

**SRM 967 – two levels**

**CAP LN 24 materials**

**SRM 909b – two levels**

### **Methods**

**NIST LC/MS method**

**U. Minnesota (J. Eckfeldt) field methods  
(CX3, Roche Jaffe, Roche Enzymatic, Vitros**

**Other Participants?**